



KaVo offers the appropriate MULTiflex®-coupling for all turbine tubing connections.

Ma.-Nr. 1.004.1387 06/02-11 en We reserve the right to make technical modifications. Slight colour differences are due to the printing process. © Copyright KaVo Dental GmbH.

Just a "Click" does the trick!



KaVo. Dental Excellence.

KaVo Dental GmbH · D 88400 Biberach/Riß · Telefon +49 7351 56-0 · Fax +49 7351 56-1488 · www.kavo.com



KaVo. Dental Excellence.

The KaVo MULTIflex®-System

Just click on - and it's ready.

Systems that are perfectly precise are both easy to use and safe. The KaVo MULTIflex®-system ensures that your work is uninterrupted when you change instruments. Regardless of the type of operation you are carrying out, the MULTIflex®-coupling will always provide you with the appropriate connection for a multitude of instruments, and at the same time it guarantees that you quickly have the right instrument at hand.

Just a quick push and a "click" – and that is sufficient to connect your new instrument onto the MULTIflex®-coupling. That means smooth operational sequences through quick changes of the instruments. What's more: the superb adaptability of the MULTIflex®-System guarantees you high efficiency through an extensive range of instruments.

When connected, the KaVo-instruments or turbines can be swivelled through a full 360°, without any annoying tension in your wrist, and without tubings becoming twisted. The instrument is simply separated from the MULTIflex®-coupling and replaced by an alternative one.

Furthermore, the built-in original KaVo non return valve prevents contaminated water from entering the system.

It's in your hands: the versatility of the MULTIflex®-coupling system saves time and money, without unnecessary interruptions of your work.

A multi functional connection.

PROPHYflex 3
polishing instrument
for professional tooth
cleaning

INTRAflex LUX
Motor instruments

Air motor 181 M and
therefore all **INTRAmatic**
shanks and **straight**
handpieces

SONICflex Airscaler,
the oscillating instrument
with a wide range of indi-
cations with SONICSYS,
SONICflex® rootplaner etc.

KaVo-turbines with
MULTIflex®-connection

KaVo RONDOflex is
the convenient alternative
to conventional
air abrasion systems.

CORONAflex for the
painless and simple
removal of crowns and
bridges via com-
pressed air

The **MULTIflex®-coupling**
is more adaptable than
you might think:

- ▶ all turbines with MULTIflex®-connection
- ▶ SONICflex, the oscillating instrument for tip operations with SONICSYS, SONICflex® rootplaner etc.
- ▶ PROPHYflex for professional tooth cleaning
- ▶ CORONAflex for the painless removal of crowns and bridges via compressed air
- ▶ all INTRAflex LUX-motor instruments
- ▶ air motor 181 M and therefore all INTRAmatic shanks and handpieces

135°C
555

All instruments can
be sterilised at up to
135° C in an autoclave.